	1 mmol/l Erdalkali-Ionen Alkaline earth ions Ions alcalino-terreux Ioni alcalino terrosi Iones alcalinos térreos Ions alcalino-terrosos	10 mg CaO/l Grad, deutsch Degrees, german Degrés, allemand Gradi tedeschi Grados, alemanes Grau alemão de dureza	10 mg CaCO <sub>3</sub> /l Grad, französisch Degrees, french Degrés, français Gradi francesi Grados, franceses Grau francês de dureza	1 mg CaCO <sub>3</sub> /l = ppm CaCO <sub>3</sub>
mmol/l	1	5,6	10	100
°dH	0,18	1	1,78	17,8
°f (TH)	0,1	0,56	1	10
ppm CaCO <sub>3</sub>	0,01	0,056	0,1	1

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
## Reagent to determine the Total Hardness in water

The Duroval offers a simple and rapid method to determine Total Hardness, based on the principle of colour change.

Each kit contains:

- 1 graduated test tube
- 1 stopper to close the test tube
- 1 tube indicator powder with a dosing spoon
- 1 bottle containing titration solution
- 1 calibrated test pipette

The indicator disperses in hard water to bright red. By adding the titration solution, the total hardness is combined, and the colour turns green. The amount of titration solution used measures the hardness contents of the water sample. The hardness value can be read off on the calibrated test pipette.

### Directions:

1. Rinse out the test tube with the water sample and fill up to the 5 ml mark (fig. 1).
2. Add 1 level dosing spoon of indicator powder and disperse by swirling (fig. 2). If the sample turns green, no hardness is present.
3. If the sample turns red, close the test tube with the stopper (fig. 3).
4. Aspirate titration solution from the storage bottle with the calibrated pipette (fig. 4 a), until the lower rim of the black piston ring is at the same level as the 0-mark (uppermost mark on the calibrated scale) (fig. 4 b). For an accurate test result, eliminate any air bubbles by depressing the plunger fully with its contents into the storage bottle, keeping the tip of the pipette immersed in the solution. Repeat the process if necessary to evacuate all air bubbles. Finally, refill the pipette to the calibrated 0-mark with the titration solution.
5. With a slight twist firmly insert the pipette tip into the orifice of the test tube stopper (fig. 5 a). Whilst shaking (fig. 5 b), depress the pipette plunger very slowly (fig. 5 c) until the colour of the water sample turns from red to green (fig. 6 a).
6. Read the total hardness taken at the lower rim of the black piston ring on the calibrated pipette (fig. 6 b) and take note. If the contents of one pipette are insufficient to effect the colour change, refill the pipette according to 4 with fresh titration solution up to the 0-mark as described and continue the titration in the same water sample according to 5 and 6.  
Note! When determining the hardness, allow for the quantity of titration solution previously used.

### Remarks:

For strongly alkaline water samples, especially boiler water e. g. pH > 10, add 1 drop of special buffer solution (order no. 400016) prior to a test (not included).

Duroval test kit Type B is suitable for the determination of Total Hardness in very soft water, i. e. below 20 mg/l (ppm).

Greater accuracy can be achieved by filling the test tube to the 10 ml mark. As twice the amount of water sample is used, the reading of the total hardness value must be divided by 2, e. g. 10 mg/l (ppm) = 5 mg/l (ppm) CaCO<sub>3</sub>.

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