# DUROVAL®

## CHLORID-TEST (Titrationsbesteck)

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## CHLORIDE TEST (Titrimetric)

With this test kit an exact determination of the chloride contents of a water sample can be carried out quickly and easily. This titration kit is an alternative to the TESTOVAL colour comparison kit, allowing determinations above 100 ppm without dilution.

Each kit contains:

- 1 dropping bottle with 18 ml reagent A
- 1 dropping bottle with 18 ml reagent B
- 2 bottles with 50 ml titration solution C
- 1 titration syringe graduated 0 - 30
- 1 test tube graduated at 5 and 10 ml
- 1 special stopper with bore hole to fit the syringe.

**Determination:** Fill the test tube exactly up to the 5 ml mark with water sample. Add 2 drops of reagent A, this will normally colour the sample blue-violet, then add reagent B dropwise until the colour changes to yellow. Close the test tube with the special stopper.

Aspire titration solution C with the syringe up to the "0"-mark. There must be no air bubble in the liquid, if an air bubble appears, it can be removed by squirting the solution back into the supply bottle by forcefully pressing down the piston, taking care that the tip of the syringe is below the level of the liquid in the bottle, and then aspiring again up to the "0"-mark.

The filled syringe is firmly set into the bore of the stopper on the test tube turning it slightly. Shaking the whole assembly gently, the piston of the syringe is slowly pressed down until the colour of the test solution just changes from yellow to violet.

The position of the lower black piston ring is read against the scale of the syringe and this value is multiplied by ten to give the correct concentration of chloride of the water sample in mg/L = ppm.

**Note:** Reagent B contains nitric acid and is corrosive. Titration solution C contains nitric acid and mercury nitrate, it is very poisonous and corrosive. Any contact with the skin should be avoided, if it should occur it must immediately be washed with plenty of tap water and soap. If the eyes or mucous membranes are affected, they should be washed with plenty of tap water as first aid, then a doctor should be consulted as fast as possible.

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