
SULFAT-TEST SULPHATE-TEST TEST du SULFATE

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With this test a quick, simple and accurate determination of sulphate in water samples can be carried out without further implements.

To achieve exact results, all metal ions must be removed from the water sample since they tend to cause errors.

For this purpose, the kit contains a small ion exchanger, which binds the metal ions and exchanges them for innocuous H⁺-ions. The resin in this exchanger was treated with a red indicator dye. When the exchanger is loaded with metal ions, a red colour gradually appears and extends slowly downwards. When after about 100 analyses 2/3 of the resin appear red, the exchanger must be regenerated.

In order to regenerate the exchanger, 20 cc of pure 10% hydrochloric acid are slowly passed in small portions through the exchanger. After the acid has been drained completely, the surplus is washed away with plenty of deionized water until no acid is found in the effluent (on testing with pH-paper, it must show over pH 4). After this regeneration, the exchanger can be used again for a similar number of analyses.

Note: After use of the exchanger, it should be flushed with 50 cc of distilled water and then closed again with cap and stopper to avoid its drying out.

Each kit contains:

- 1 Ion exchanger with cap and stopper
- 2 bottles with 50 cc each of reagent "A"
- 1 dropper bottle with 8 cc of reagent "B"
- 1 bottle with 50 cc of titrant "C"
- 1 graduated glass titration vial with connective stopper and perforated screw cap
- 1 titration syringe
- 1 test tube, graduated at 1 and 10 cc for dilution.

Procedure:

- 1.) Cap and stopper are removed from the exchanger. About 10 cc of the water sample are run through and discarded (fig. 1).
- 2.) Another portion of the same water is poured through and of this effluent 5 cc are collected in the glass vial until it is filled exactly up to the lowest mark (fig. 2).
- 3.) 4 cc of reagent A are added until the next mark is reached (fig. 3).
- 4.) One drop of reagent B is added (fig. 4).
- 5.) The vial is closed with the connective stopper and the latter secured by screwing on the perforated cap (fig. 5).
- 6.) The syringe is filled with titrant C up to the "0" mark (fig. 6). If air bubbles appear in the cylinder, the piston is pressed forcefully all the way down and then the liquid aspirated again.
- 7.) The syringe is pressed securely into the bore of the connective stopper (fig. 7) and with gentle shaking the piston is very slowly pressed downwards. The first minimal amount of titrant will already effect a colour change from violet to light blue. If sulphate is present, the colour will after a few seconds revert to violet. The titration is continued in this manner until the light blue colour persists for at least one minute.
- 8.) When the blue colour persists, the value is read from the position of the black piston ring on the graduation of the cylinder of the syringe (fig. 8). The value read there, multiplied by ten = ppm sulphate (e.g. 8 = 80 mg SO₄⁻²/L). If the contents of one syringe are not sufficient to effect the permanent colour change, (300 ppm SO₄⁻²), the syringe can be filled again with titrant up to the "0" mark and the titration continued in the same manner until the light blue colour persists. The value found then is added to the 300 ppm of the first syringe.

If even the contents of the second syringe are not sufficient to effect a permanent colour change, the water sample should be diluted 1 : 10 with deionized water and the entire analyse repeated completely, including the passage of the exchanger with this diluted sample.

The value read then must be multiplied by 100.

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