testoval®

AMMONIUM-TEST

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Directions for Use

Remove the test cell from the comparator and fill it up to the mark with water sample. Add 5 drops of reagent A and dissolve by stirring with the red marked spoon. Then add 1 spoon of reagent B and stir again. After 5 minutes add 10 drops of reagent C and stir, 5 minutes later, reinsert the test cell into the comparator, hold the comparator against a bright background or light source, and match the colour intensity of the central field (test solution) with one of the colours in the peripheral colour fields. Read the value in ppm NH_4^+ , which is printed under the matching colour field.

Should the colour in the central field be more intensive than any of the colours in the peripheral fields, repeat the test with a water sample diluted 1 : 10. For this purpose fill the graduated plastic tube up to the lower mark with sample water and dilute with ammonium-free water (e.g. distilled water) to the upper mark. Transfer the diluted sample into the test cell, which has been rinsed with ammonium-free water and proceed as described above. The value obtained must in this case be multiplied by 10.

Caution: In case of ammonium-free water colour changes to intensive yellow. The reaction must be carried out at room temperature (20°C - 25°C).

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